

New!

from *Mada*

No. 1702D



No. 1702B



No. 1702 Series Flowmeters

No. 1702D



No. 1702B



The Mada Flowmeter is designed for use on hospital and emergency oxygen equipment where a regulated supply of oxygen is required. The Mada Flowmeter has 11 flow settings that are highly accurate and easy to read.

The durable, low profile design allows this product to be used under rugged conditions. The inlet thread is 1/8" NPT to allow the use of standard oxygen connections.

Available Models:

No. 1702D Oxygen Flowmeter with standard male DISS outlet.

No. 1702B Oxygen Flowmeter with standard barb outlet.

Oxygen Inlet Adapters:

No. 1703 DISS Female. **No. 1706** Ohio Quick Connect.



MADA, INC.

DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

See our Web Site: www.madamedical.com



Mada Flowmeter

Designed for use on hospital and emergency oxygen equipment where a regulated supply of oxygen is required.

The Mada Flowmeter has 11 flow settings that are highly accurate and easy to read. Flow rates: .5, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25 LPM.

Easy to turn liter flow knob that locks into place for each liter setting. Liter flow numbers are large enough to see in low light conditions.

The durable, low profile design allows this product to withstand every day use.

The inlet thread is 1/8" NPT allowing the use of standard oxygen connections.

Aluminum and polycarbonate construction provides a Flowmeter that is both lightweight and durable.

Available in several configurations to handle most applications.

Specifications

Inlet Pressure: 50 psig nominal

Operating Temperature: -30° F to 125° F

Storage Temperature: -40° F to 160° F

Flow Rates: .5, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25 LPM

Filter: 40 micron porous bronze

Inlet Port: 1/8" NPT female

Outlet Fittings: Standard Hose Barb
Male DISS fitting

Accuracy: Plus or minus 10%

Material: Body-Anodized aluminum
Knob-Polycarbonate



MADA, INC.

DOMESTIC/INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072

Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

Web Site: www.madamedical.com

Distributed By: