



New!

from *Mada*

Multi-O2-Manifold

No. 1713M

The Mada Multi-O2-Manifold is a compact O2 Administrator with 8 flow control valves mounted on a collapsible tripod. The Multi-O2-Manifold can be connected in series with a second unit that can expand the capabilities to treating 16 patients..

- No assembly necessary.
- Selector valves provide .5 - 25 LPM.
- Allows for administration of low flow oxygen and/or oxygen powered resuscitation.
- Each unit comes equipped with an extra self closing port for a demand valve or aspirator.

Available Models:

No. 1713M Multi-Flow-O2-Manifold.

No. 1714M Multi-Flow-O2-Manifold with Humidifiers, Masks and 10' Hose.

Optional Accessories

1534MA Demand Valve with Adult Mask and 6' Hose.

1543 Aspirator with 6' Hose and Catheter.

1327 Adult size Mask and 7' Tube (Box/50).

1543-H-10 10' Conductive Hose.

302-4P Plastic Hose Connector.

Specifications

Supply Pressure: 50 psig nominal

Operating Temperature: -30°F to 125°F

Storage Temperature: -40°F to 160°F

Filter: 40 Micron

Material: Body: Anodized Aluminum
Fittings: Plated Brass

Outlets: 9 outlets including 8 Flow Controllers with flow rates of .5, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25 LPM and one standard male DISS check valve outlet

Distributed By:



MADA, INC.

DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509
To Order Call Toll Free: **800-526-6370**
Web Page: www.madamedical.com

MFOM 6/05