

New!

from *Mada*

Disposable BVM Resuscitator



- Medication port for use with MDI or syringe.
- Manometer port.
- Integrated PEEP connection, no adapter necessary.
- Adult, Child and Infant sizes.
- Pressure release valve standard on child and infant models, optional on adult model.
- All models complete with reservoir bag and O₂ tubing.
- Transparent bag (green) with positive textured grip and sensitive response.
- Range of accessories including 3M™ HEPA Filter and PEEP valve.
- Latex free.
- Disposable.



MADA, INC.

DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

Web Page: www.madamedical.com