

No. 7000-MVP

Mada

Disposable Manual Vacuum Pump



- *State of the art, high vacuum, light, portable emergency suction device.*
- *Built-in overflow protection seal.*
- *100° swivel allows better accessibility in tight situations.*
- *Efficient operation and simple maintenance.*
- *One hand operation.*
- *Large volume capacity.*
- *Disposable collection jar and catheter.*



MADA Medical Products, Inc.
DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call: **800-526-6370**

Web Page: www.madamedical.com

r 6/10



Technical Information

<i>Vacuum (Max):</i>	-550mmHg
<i>Peak free airflow:</i>	30 l/min
<i>Disposable container volume:</i>	300ml
<i>External diameter of cap (Suction Port):</i>	17mm
<i>External diameter of catheter connection:</i>	13mm and 6mm
<i>Operating environmental temperature:</i>	-20°C + 50°C
<i>Storage environmental temperature:</i>	-40°C + 60°C

*Performance values given are achievable under test conditions.
Values may vary during actual use.*

Ordering Information

No. 7000-MVP *Mada Manual Vacuum Pump, Disposable Collection Jar with adult and child catheter*

No. 7000-AD *Disposable Collection Jar with adult and child catheter*



MADA Medical Products, Inc.

DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call: **800-526-6370**

Web Page: www.madamedical.com