

MADA DISPOSABLE OXYGEN MASKS AND ACCESSORIES

All Products Latex Free

Mada emergency oxygen masks and accessories are made from high quality materials which result in a soft, clear finish for greater patient comfort and easy monitoring of patient condition.



No. 1327 - Adult with 7' Oxygen Tubing, Medium Concentration, case of 50.



No. 1327C - Child with 7' Oxygen Tubing, Medium Concentration, case of 50.



No. 5104 - Adult with 7' Oxygen Tubing, High Concentration, Rebreathing with Reservoir, case of 50.



No. 53224 - Child with 7' Oxygen Tubing, High Concentration, Rebreathing with Reservoir, case of 50.



No. 5108 - Adult Non-Rebreathing Mask with 7' Oxygen Tubing, with Reservoir and two Side Valves, case of 50.



No. 1326 - Adult Partial Non-Rebreathing Mask with 7' Oxygen Tubing, with Reservoir and one Side Valve, case of 50.



No. 53226 - Child Partial Non-Rebreathing Mask with 7' Oxygen Tubing, with Reservoir and one Side Valve, case of 50.



No. 53266 - Infant, High Concentration, Non-Rebreather with 7' Oxygen Tubing, Check Valve, two Side Valves, case of 50.



No. 1325-7 - Adult Nasal Cannula with Curved Tip, 7' "No-Crush" Oxygen Tubing, case of 50.