

# MADA EMERGENCY OXYGEN KITS

No. 1301ME



No. 1301BME



No. 1527E



No. 1528E



No. 1530E



No. 1528BE



No. 1531E



No. 1533E



No. 1631E



No. 1629E



**MADA, INC.**

DOMESTIC/ INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072

Toll Free: 800-526-6370 Tel: 201-460-0454 Fax: 201-460-3509

See our Web Page: [www.madamedical.com](http://www.madamedical.com) for a listing of our products.

To order unit filled, remove suffix "E"

EO<sub>2</sub> KITS 6/05

# MADA EMERGENCY OXYGEN UNITS

<b>PRODUCT NO.</b>	<b>CYLINDER</b>	<b>REGULATOR</b>	<b>CASE/CART</b>	<b>FEATURES</b>
1201ME	(M7) 1203ME	1308A 6 LPM Avg.	1218A	Inhalator only
1301ME	(M7) 1303ME	1308A 6 LPM Avg.	1318A	Inhalator only
1300ME	(M7) 1303ME	1308A 6 LPM Avg.	1319A	Inhalator with extra cylinder
1301BME	(M7) 1303ME	1308A 6 LPM Avg.	1319A	With BMR
1331ME	(M7) 1303ME	1308A 6 LPM Avg.	1319A	With BMR and extra cylinder
1501E	(MD) 1503E	1442 6 LPM Avg.	1520	Inhalator only
1501ME	(MD) 1503ME	1308A 6 LPM Avg.	1520	Inhalator only
1527E	(MD) 1503ME	1308A 6 LPM Avg.	1527C	Wall Case Inhalator
1528E	(MD) 1502E	1441 6 LPM Avg.	1500A	Inhalator only
1528BE	(MD) 1502E	1441 6 LPM Avg.	1500A	Inhalator with BMR
1529E	(MD) 1502E	R1835-25GB (1/4-25 LPM)	1520	Inhalator only
1530E	(MD) 1502E	R1835GBHF-25 (1/4-25 LPM)	1500A	Demand Valve, with Inhalator
1531E	(MD) 1502E	1635-1 (50 PSI Reducer)	1500A	Demand Valve, with no Inhalator
1533E	(MD) 1502E	R1835GBHF-25 (1/4-25 LPM)	1500A	Deluxe with aspirator
1629E	(ME) 1603E	1442 6 LPM Avg.	370C	Inhalator with BMR
1629DE	(ME) 1603E	1442-10 10 LPM Avg.	370C	Adult/Child BMR
1630ADE	(ME) 1602E	1335 (2-8 LPM)	370C	Disposable BMR
1631E	(ME) 1602E	1635-1 (50 PSI Reducer)	370C	Demand Valve, no Inhalator
1634E	(ME) 1602E	R1835GBHF-25 (1/4-25 LPM)	370C	Demand Valve, with Inhalator

Other valve combinations available upon request.

Filled Kits available upon request.

Note: BMR is abbreviation for Bag Mask Resuscitator used for manual resuscitation.

Demand Valve is used for Oxygen powered Resuscitation.

Inhalator is for administering low flow oxygen only.



**MADA, INC.**

DOMESTIC/ INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072

Toll Free: 800-526-6370 Tel: 201-460-0454 Fax: 201-460-3509

See our Web Page:

[www.madamedical.com](http://www.madamedical.com) for a listing of our products.

To order unit filled, remove suffix "E"

Distributed By:

EO<sub>2</sub>KITS 6/05