


NO. 1411E

*Mada*

# First-In Oxygen Kit



*The New Mada First-In Oxygen Kit is designed to be used by first responders, EMS agencies, fire departments, law enforcement agencies and occupational health organizations.*

*The kit contains enough equipment to care for a patient in respiratory distress or respiratory arrest.*

*With the optional AED, it can be used to begin treatment of a patient suffering from sudden cardiac arrest.*

*The kit is designed so that the care giver can begin to treat the patient while waiting for an ambulance to arrive.*

*Each kit comes complete with a 248 liter Luxfer aluminum oxygen cylinder (Empty) with built-in on/off toggle valve, oxygen regulator (1/2-15 LPM), disposable adult bag mask resuscitator, adult nasal cannula, adult non-rebreather mask and a padded carrying bag with padded shoulder strap and easy carry handle.*

## **Ordering Information:**

**No. 1411E:** First-In Oxygen Kit, includes aluminum cylinder, (Empty) oxygen regulator (1/2-15 LPM), nasal cannula, non-rebreathing mask, disposable adult bag mask resuscitator and carrying bag.

**Note:** To have Oxygen Kit filled, remove suffix "E".



**MADA, INC.**

DOMESTIC/INTERNATIONAL  
625 Washington Avenue, Carlstadt, NJ 07072  
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

See our Web Site: [www.madamedical.com](http://www.madamedical.com)



## Versatile

Kit contains enough equipment to care for a patient in respiratory distress or respiratory arrest. Care giver can begin to treat the patient while waiting for the ambulance to arrive.

## Mada Quality

All components are manufactured to the exceptional quality Mada is known for.

## Durable

Product is designed to hold up under EXTREME conditions.

## Complete Kit

Includes everything needed to provide oxygen and ventilation.

## Specifications

**Supply Pressure:** 500-2,300 PSI

**Flow Rate of Regulator:** .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM.

**Bag Mask Resuscitator:** Disposable, Single Patient Use.  
Adult, with O<sub>2</sub> reservoir and mask.

**Weight:** Approximately 8 pounds.

## First-In Kit No. 1411E Includes:

**No. 1402E** Aluminum O<sub>2</sub> cylinder. {Empty}

**No. M1935-15GB** O<sub>2</sub> regulator (1/2 - 15 LPM).

**No. 4000X** Disposable Adult Bag Mask Resuscitator.

**No. 1325-7X** Nasal Cannula with 7' tube.

**No. 1326X** Non-rebreather Mask and 7' tube.

**No. 1410** Padded Carry Bag with padded strap.



### MADA, INC.

DOMESTIC/INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072

Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

See our Web Site: [www.madamedical.com](http://www.madamedical.com)

Distributed By: