

MADACYLINDER™ III

EMERGENCY & THERAPY

OXYGEN UNITS

Featuring the MadaValve System



No. 1303ME

Over thirty years ago Mada Medical introduced its MadaCylinder™ line.

Today, a third generation, the MadaCylinder™ III offers a new dimension in design.

Exclusively from Mada, the MadaCylinder™ III is ultra light, weighing only 4 1/4 pounds and contains approximately 209 liters of oxygen.

The MadaCylinder™ III is available in a variety of portable oxygen units to meet the needs of dealers and their patients.

Features:

- Ultra-light-weight, cylinder weight 4 1/4 lbs.
- Volume of oxygen, 209 liters.
- Flat bottom design.
- Heat sensitive coating.
- Built-in contents gauge. (Mada Innovation)
- Refillable.
- Built-in pressure and temperature safety devices.
- Uses MadaValve regulators.
- Many models to choose from to meet your patients' needs.

Benefits:

- Small and easy to carry.
- Long lasting, up to approximately 100 minutes at 2 LPM.
- Indicates volume in cylinder even when not in use.
- Uses standard CGA-540 outlet thread.
- Safe to use, meets or exceeds all government specifications.
- Compatible with Mada regulators No. 1308A and No. 1309A series.
- Available in single cylinders or kits.



No. 1315ME

1315ME Oxy-Uni-Pak

Contains: one 1303ME MadaCylinder III, 1309A adjustable flow regulator (2 - 8 LPM), 7 ft. nasal cannula, in shoulder bag.

1315M-5E - Same as 1315ME except uses 1309A-5 regulator (1/2 - 5 LPM).

Therapy Oxygen Units



No. 1314 ME

1314ME Oxy-Uni-Pak

Contains: one 1303ME MadaCylinder III, 1309A adjustable flow regulator (2 - 8 LPM), 7 ft. nasal cannula, in carry case.



No. 1313ME

1313ME Oxy-Duo-Pak

Contains: two 1303ME MadaCylinder III's, 1309A adjustable flow regulator (2 - 8 LPM), 7 ft. nasal cannula, in carry case.

To order units filled, remove suffix "E."

EO₂TOU 6/05

EMERGENCY UNITS



No. 1301ME Oxy-Uni-Pak

A complete single cylinder oxygen unit in custom carrying case.
Contains: one 1303ME MadaCylinder III, 1308A fixed flow (6 LPM Avg.) regulator, mask and tube, in carry case.



No. 1300ME Oxy-Duo Pak

A complete emergency oxygen unit.
Contains: two 1303ME MadaCylinder III, 1308A fixed flow (6 LPM Avg.) regulator, mask and tube, in carry case.



No. 1301BME C.P.R. Pak

A compact, emergency oxygen resuscitation unit.
Contains: one 1303ME MadaCylinder III, 1308A fixed flow (6 LPM Avg.) regulator, 1432 resuscitator, mask and tube, in carry case.



No. 1303ME MadaCylinder III

Ultra-lightweight oxygen cylinder.
Built in contents gauge.



No. 1331ME C.P.R. Pak

A lightweight portable unit with manual resuscitator.
Contains: two 1303ME MadaCylinder III, 1308A fixed flow (6 LPM Avg.) regulator, mask and tube, carry case, and 1432 manual resuscitator, packed separately.



PLUS ►



No. 1308A

Fixed flow regulator (6 LPM Avg.) with CGA-540 thread, this regulator screws on easily by hand to MadaValve Cylinders.



No. 1309A

Adjustable flow regulator (2-8 LPM) with CGA-540 thread, screws on easily by hand to MadaValve Cylinders.
Also available as **1309A-1** 1/16-1 LPM, **1309A-5** 1/2-5 LPM and **1309A-15** 2-15 LPM.



MADA, INC.

DOMESTIC/ INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Toll Free: 800-526-6370 Tel: 201-460-0454
Fax: 201-460-3509

Web Page: www.madamedical.com

To order unit filled, remove suffix "E"

Distributed By: