


MadaVac 172BS-II Aspirator

from **MADA!**

MadaVac 172BS-II
Heavy duty Aspirator...Designed with a
high performance vacuum pump. Created for
home care and institutional use.



Standard Features

- Bacteria filter prevents contaminant damage.
- Vibration free, easy to read gauge.
- Quick connect tubing for all connections.
- Draws up to 22" Hg (560mmHg) with a flow rate up to 40 LPM.
- Large 800cc collection bottle with overflow trap.
- Ergonomically designed carrying handle.
- Quiet, yet powerful 115V/60Hz/0.68A
- Two year warranty.

Common Uses: Hospitals, Nursing Facilities, Physician Offices, Homecare.

Accessories

- 178B** - 800cc Suction Bottle, single use with overflow trap (*Box 10*).
- 178TB** - 6' Suction Tubing (*Box 10*).
- 178BE** - Elbow Fitting (*Each*).
- 3030200040** - 7" Vacuum Bottle's Silicone Tube (from filter to bottle) (*Each*).
- 3030200120** - 5.5" Silicone Tube (from filter to regulator) (*Each*).
- A05020030** - Bacteria Filter (*Each*).
- 172VGTEST** - Vacuum test gauge assembly (certified) (*Each*).



MADA, INC.
DOMESTIC/INTERNATIONAL
625 Washington Ave., Carlstadt, NJ 07072
Tel: 201-460-0454 Toll Free: 800-526-6370
Fax: 201-460-3509