

No. 7000-MVP

**New!**

from *Mada*

# Disposable Manual Vacuum Pump



- *State of the art, high vacuum, light, portable emergency suction device.*
- *Built-in overflow protection seal.*
- *100° swivel allows better accessibility in tight situations.*
- *Efficient operation and simple maintenance.*
- *One hand operation.*
- *Large volume capacity.*
- *Disposable collection jar and catheter.*



**MADA, INC.**

DOMESTIC/INTERNATIONAL  
625 Washington Avenue, Carlstadt, NJ 07072  
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

Web Page: [www.madamedical.com](http://www.madamedical.com)



## Technical Information

<i>Vacuum (Max):</i>	<i>-550mmHg</i>
<i>Peak free airflow:</i>	<i>30 l/min</i>
<i>Disposable container volume:</i>	<i>300ml</i>
<i>External diameter of cap (Suction Port):</i>	<i>17mm</i>
<i>External diameter of catheter connection:</i>	<i>13mm and 6mm</i>
<i>Operating environmental temperature:</i>	<i>-20°C + 50°C</i>
<i>Storage environmental temperature:</i>	<i>-40°C + 60°C</i>

*Performance values given are achievable under test conditions.  
Values may vary during actual use.*

## Ordering Information

**No. 7000-MVP** Mada Manual Vacuum Pump, Disposable Collection Jar with adult and child catheter

**No. 7000-AD** Disposable Collection Jar with adult and child catheter



### **MADA, INC.**

DOMESTIC/INTERNATIONAL  
625 Washington Avenue, Carlstadt, NJ 07072  
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

Web Page: [www.madamedical.com](http://www.madamedical.com)