



Mini-Oxygen-Manifold

No. 1712M

The Mada Mini-Oxygen-Manifold is an easy to operate constant flow device that provides oxygen for up to five patients at the same time.

Several Mada Mini-Oxygen-Manifolds can be connected for emergency situations where more oxygen outlets are required.

Available Models:

No. 1712M Mini-Oxygen-Manifold with 5 Oxygen DISS check valves, (5) 15 LPM orifices and carry bag.

No. 1712M-6 Mini-Oxygen-Manifold with 5 Oxygen DISS check valves, (5) 6 LPM orifices and carry bag.

No. 1712M-10 Mini-Oxygen-Manifold with 5 Oxygen DISS check valves, (5) 10 LPM orifices and carry bag.

Five Oxygen Outlets

These outlets can supply 40-60 psig or be used with the low flow orifices supplied.

Versatile

Can provide 50 psig to drive demand valves, resuscitators, flow meters, or suction devices.

Expandable

Can be expanded for emergency situations where more oxygen outlets are required.

Simple Design

Offers standard DISS outlets. Easy to connect orifices.

Specifications

Supply Pressure: 40-60 psig

Operating Temperature: -30°F to 125°F

Storage Temperature: -40°F to 160°F

Flow Rates: 6, 10, or 15 LPM

Weight: 13.5 oz.

Sizes: Without Barbs: 1" x 3" x 6"

With Barbs Attached: 1" x 6" x 7 1/2"

Material: Anodized aluminum

Manifold: Anodized aluminum

Fittings: Plated brass

Barbs: H. I. polystyrene

Optional Orifices

No. 302-4PW-6 6 LPM orifice pkg. 5

No. 302-4PW-10 10 LPM orifice pkg. 5

No. 302-4PW-15 15 LPM orifice pkg. 5 (included with 1712M)



MADA, INC.

DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072

Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

See our Web Site: www.madamedical.com