

# Mada Oxygen Pressure Regulators



*Mada...  
Quality  
You Can  
Count On!*

*Mada offers a wide range of oxygen regulators that are designed for both hospital and EMS applications at the most competitive prices.*

*Mada oxygen regulators offer the features you expect, including easy-to-read flow rates and a dial-in click-style flow controller.*

**The Mada series of all-brass body oxygen pressure regulators provide superior safety features and a greater number of flow selections for the EMS professionals.**

## **Features**

- *100% Brass construction for durability and safety.*
- *Large, easy-to-read flow rate numbers.*
- *Compact, low-profile contents gauge with gauge protector.*



**MADA, INC.**

DOMESTIC/INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072  
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

See our Web Site: [www.madamedical.com](http://www.madamedical.com)

# Specifications

---

<b>Supply Pressure:</b>	200-2600 psi
<b>Outlet Pressure:</b>	High flow outlets: 45-55 psi at 200-2,600 psi inlet pressure
<b>Flow Capacity:</b>	High flow outlets: 100 lpm minimum single or combined flow with 200 psi inlet pressure Constant flow outlet selections: 0-25 lpm: 0, 1/4, 1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25
<b>Filter:</b>	25 micron sintered bronze
<b>Relief Valve:</b>	Yes
<b>Operating Temperature:</b>	-20° F to 125° F
<b>Storage Temperature:</b>	-20° F to 125° F
<b>Materials:</b>	Body: Nickel-plated brass Outlets: Nickel-plated brass Internal parts: Brass, stainless steel, Viton®, silicon and Teflon®

## Ordering Information

---

### **Yoke Pin Indexed CGA-870 Model Regulators For Small Oxygen Cylinders.**

<b>R1835GBHF-25</b>	1/4-25 lpm with gauge, 2 DISS check valves for 50 psi hose barb outlets.
<b>R1835-25BR</b>	1/2-25 lpm with gauge, hose barb outlet.
<b>1635-1</b>	Pressure reducer with gauge, 1 DISS 50 psi outlet.

### **CGA-540 Model Regulator For Large Oxygen Cylinders.**

<b>R1833GBHF-25</b>	1/4-25 lpm with gauge, 2 DISS check valves for 50 psi hose barb outlets.
<b>1635-1</b>	Pressure reducer with gauge, 1 DISS 50 psi outlet.

**Additional models are available. For information about our complete line of oxygen regulators. Please contact your Mada customer service representative or visit us online.**



**MADA, INC.**

DOMESTIC/INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072  
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

See our Web Site: [www.madamedical.com](http://www.madamedical.com)

10/11