

*New*

# Penguin Compressor/Nebulizer

from **MADA!**



## Standard Features

- Powerful piston pump.
- Lightweight: 4 lbs. 2 oz.
- Child friendly.
- Built in handle.
- Maximum pressure: 35 PSI.
- Operating pressure: 18 PSI.
- Liter flow: 6-7 LPM at 11 PSIG.
- Comes complete with nebulizer, mouthpiece, 7' tubing, and mask.
- Padded shoulder bag included.
- 3 year warranty.
- Electrical requirements: 115 volts, 60Hz.

## Accessories

- 5752** - Replacement Medication Nebulizer with tee adapter, mouthpiece and 7' tubing. (*Box 50*).
- 5754** - VixOne Nebulizer, Pacifier, 45 degree elbow and 7' tubing. (*Box 50*).
- 50311** - Replacement Medication Nebulizer with pediatric mask and 7' tubing. (*Box 50*).
- 277AF** - Replacement Air Filters (*Box 5*).



**MADA, INC.**

DOMESTIC/INTERNATIONAL  
625 Washington Avenue., Carlstadt, NJ 07072  
Tel: 201-460-0454 Toll Free: 800-526-6370 Fax: 201-460-3509

See our web page: [www.madamedical.com](http://www.madamedical.com)