


Penguin Compressor/Nebulizer

from **MADA!**

No. 278

Penguin

278SB *Optional:*
Padded houlder bag



The Penguin is a child friendly Compressor/Nebulizer that provides quieter treatment and durability

Standard Features

- Powerful piston pump.
- Lightweight: 4 lbs. 2 oz.
- Child friendly.
- Built in handle.
- Maximum pressure: 35 PSI.
- Operating pressure: 18 PSI.
- Liter flow: 6-7 LPM at 11 PSIG.
- Comes complete with nebulizer, mouthpiece, 7' tubing, and mask.
- 3 year warranty.
- Electrical requirements: 115 volts, 60Hz.

Accessories

- 5752 - Replacement Medication Nebulizer with tee adapter, mouthpiece and 7' tubing. (Box 50).
- 5754 - VixOne Nebulizer, Pacifier, 45 degree elbow and 7' tubing. (Box 50).
- 50311 - Replacement Medication Nebulizer with pediatric mask and 7' tubing. (Box 50).
- 277AF - Replacement Air Filters. (Box 5).
- 278SB - Padded shoulder bag with strap.



MADA, INC.

DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509
To Order Call Toll Free: **800-526-6370**
Web Site: www.madamedical.com

Quality products since 1969