

Mada

No. 1432



Silicone BVM Resuscitator



- Completely autoclavable and reusable.
- Textured silicone bag with positive grip and sensitive response to compliance.
- Versatile 360° swivel for facemask or endotracheal tube.
- Child bag with 40cm safety pressure release valve.
- Clear non-rebreathing patient valve.
- Folds into compact convenient size.
- Full range of accessories (Adult, Child size facemasks, PEEP valve).
- Designed to meet A. S. T. M. and European standards.
- Durable silicone bag and mask with polycarbonate fittings.
- Complete with reservoir bag and O₂ tubing.
- Latex free.



MADA Medical Products, Inc.

DOMESTIC/INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call: **800-526-6370**

Web Page: www.madammedical.com



Technical Specifications

	<i>Infant Bag</i>	<i>Child Bag</i>	<i>Adult Bag</i>
<i>Bag volume</i>	<i>280ml</i>	<i>550ml</i>	<i>1,500ml</i>
<i>Oxygen reservoir volume</i>	<i>600ml</i>	<i>2,500ml</i>	<i>2,500ml</i>
<i>Length</i>	<i>15cm</i>	<i>15cm</i>	<i>23cm</i>
<i>Folded length</i>	<i>13cm</i>	<i>7.5cm</i>	<i>13cm</i>
<i>Weight</i>	<i>295g</i>	<i>363g</i>	<i>544g</i>

Connectors: 15mm I. D. and 22mm O. D.

Operation temperature: -18°C to 50°C

Storage temperature: -40°C to 60°C

Materials (All components are latex-free)

Bag and mask: Silicone

Valve and fittings: Polycarbonate

Lip and butterfly valves: Silicone

Oxygen reservoir and tubing: PVC

Ordering Information

No. 1428 Reusable Bag Mask Resuscitator, Infant

No. 1432 Reusable Bag Mask Resuscitator w/reservoir, Adult

No. 1434 Reusable Bag Mask Resuscitator w/reservoir, Child

No. 1435 Reusable Adult Mask for No. 1432

No. 1436 Reusable Child Mask for No. 1434

No. 1436NN Reusable Infant Mask for No. 1428

No. 1439 Case for Manual Resuscitator



MADA Medical Products, Inc.

DOMESTIC/INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072

Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call: **800-526-6370**

Web Page: www.madamedical.com