


The Mada

ULTRASONIC *Cleaning System*



Designed specifically for cleaning the MadaJet.[®]

Includes custom cover and rack that holds 2 MadaJets.

The Mada Ultrasonic Cleaning System provides super-strength cleaning every time. Each tank and body is carefully inspected for strength and durability.

FEATURES:

- Increased power, yet cleans quietly
- Perfect for countertop use
- 60 minute timer
- Includes custom lid and rack for the MadaJet
- Design fits into any Medical/Dental Laboratory
- Backed by 2 1/2 year warranty



MADA, INC.
DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509
Toll Free: 800-526-6370

The Mada **ULTRASONIC** Cleaning System

DIMENSIONS:


Overall 10.25" x 6.5" x 8.25"

Internal: 9.38" x 5.38" x 4.0"

Weight: 10 Pounds

Tank Capacity: 3.38 Quarts

CE Approved: Yes




No. 6008

Order Information:

- 6008** *Ultrasonic Unit with custom lid and rack.*
- 6009** *Rack.*
- 6010** *Universal Basket.*
- 7021** *MadaCide-FD case of 4 Gallons for disinfecting MadaJet and accessories.*



No. 7021



No. 6009

No. 6010

To Order Call Toll Free

800-526-6370

See our web site

www.madamedical.com



MADA, INC.
DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 Fax: 201-460-3509
Toll Free: 800-526-6370