

Mada's Voyage-II Compressor/Nebulizer



*Voyage-II
Delivers effective
and efficient aerosol
medication therapy
with clinical accuracy.
Built for quality and durability
for years of reliable treatments.*

No. 277

Standard Features

- Powerful piston pump.
- Lightweight: 4 lbs. 2 oz.
- Storage for accessories.
- Built in handle.
- Maximum pressure: 35 PSI.
- Operating pressure: 18 PSI.
- Liter flow: 6-7 LPM at 11 PSIG.
- Comes complete with nebulizer, tee adapter, mouthpiece and 7' tubing.
- Comes with (5) replacement air filters.
- 3 year warranty.
- Electrical requirements: 115 volts, 60Hz.

Accessories

- 5752** - Replacement Medication Nebulizer with tee adapter, mouthpiece and 7' tubing. (Box 50).
- 1310NMA** - Replacement Medication Nebulizer with adult mask and 7' tubing. (Box 50).
- 50311** - Replacement Medication Nebulizer with pediatric mask and 7' tubing. (Box 50).
- 277AF** - Replacement Air Filters (Box 5).



MADA, INC.
DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509
To Order Call Toll Free: **800-526-6370**
Web Site: www.madamedical.com

Quality products since 1969