

Mada's **Voyage-II** **Compressor/Nebulizer**



Voyage-II

*Delivers effective
and efficient aerosol
medication therapy
with clinical accuracy.
Built for quality and durability
for years of reliable treatments.*

No. 277

Standard Features

- Powerful piston pump.
- Lightweight: 4 lbs. 2 oz.
- Storage for accessories.
- Built in handle.
- Maximum pressure: 35 PSI.
- Operating pressure: 18 PSI.
- Liter flow: 6-7 LPM at 11 PSIG.
- Comes complete with nebulizer, tee adapter, mouthpiece and 7' tubing.
- Comes with (5) replacement air filters.
- 3 year warranty.
- Electrical requirements: 115 volts, 60Hz.

Accessories

- 5752** - Replacement Medication Nebulizer with tee adapter, mouthpiece and 7' tubing. (*Box 50*).
- 1310NMA** - Replacement Medication Nebulizer with adult mask and 7' tubing. (*Box 50*).
- 50311** - Replacement Medication Nebulizer with pediatric mask and 7' tubing. (*Box 50*).
- 277AF** - Replacement Air Filters (*Box 5*).



MADA, INC.

DOMESTIC/INTERNATIONAL
625 Washington Avenue., Carlstadt, NJ 07072
Tel: 201-460-0454 Toll Free: 800-526-6370 Fax: 201-460-3509

See our web page: www.madamedical.com