

IT'S ALL IN OUR NAME!

MadaCide-FD

FAST DRYING (FD) and more...

Now KILLS
SWINE FLU
(H1N1)



✓ One Step
Infection
Control Spray-
Air Dry

✓ 3 Minute MRSA
Kill Time

✓ 5 Minute
VRE Kill Time

✓ 6 Minute
TB Kill Time

✓ 30 Second
HIV Kill Time

✓ 3 Minute
Hepatitis B & C
Kill Time



MADA Medical Products, Inc.
625 Washington Avenue
Carlstadt, NJ 07072
800-526-6370
Fax: 201-460-3509

MadaCide-FD HOSPITAL DISINFECTANT/CLEANER
MEETING YOUR INFECTION CONTROL NEEDS FROM MADA
A NAME YOU HAVE TRUSTED FOR OVER 42 YEARS

MadaCide-FD

(FAST DRYING)

Another weapon in the fight for Infection Control

YOU ASKED FOR IT!

HERE'S HOW FAST IT WORKS

ORGANISMS	MINIMUM CONTACT TIME
Methicillin Resistant Staphylococcus Aureus (MRSA)	3 minutes
Mycobacterium bovis BCG (Tuberculosis)	6 minutes
Hepatitis B & C Virus (DHBV), Polio I Virus, Pseudomonas aeruginosa, Rhinovirus, Salmonella choleraesuis and Staphylococcus aureus	3 minutes
Aspergillus niger, Candida albicans, Vancomycin-resistant Enterococcus faecalis (VRE), Canine Parvovirus and Trichophyton mentagrophytes	5 minutes
Herpes Simplex II, HIV-I (Aids Virus) and Influenza Strain A2/HK	30 seconds

MadaCide-FD is EPA registered, assuring that all performance claims stated on the label are accurate and verifiable.

It is Fungicidal, Pseudomonacidal, Virucidal, Bactericidal and Tuberculocidal.

Its unique anti-corrosion formulation will not damage lenses, cement, plastics, rubber, steel, aluminum, brass or any other metal when used as directed.

For Surgical Instrument Presoak:

Remove gross debris. Presoak surgical instruments in undiluted solution for 10 minutes. Proceed with normal decontamination procedure.

AREAS OF USE:

*Hard surfaces in: Renal Dialysis Centers,
Surgery, Recovery, Anesthesia, X-Ray,
Cath. Lab, E. R., Orthopedics,
Newborn Nursery, Respiratory Therapy,
Radiology, Central Supply,
Ultrasonic Cleaners*

ORDERING INFORMATION:

No. 7020 - 32 oz. Bottle Spray..... 12/case
No. 7021 - 1 Gallon..... 4/case
No. 7022 - 2.5 Gallon..... each
No. 7023 - 5 Gallon..... each
No. 7024 - 55 Gallon..... each



MADA Medical Products, Inc.
625 Washington Avenue
Carlstadt, NJ 07072
800-526-6370
Fax: 201-460-3509

Distributed By:

See our Web Site: www.madamedical.com for a listing of our products