


MadaMeter

from *Mada*

An oxygen flow tester at your fingertips!

The MadaMeter is a portable hand held flowmeter which tests oxygen flow rates quickly and easily.

- *Easy to use.*
- *Silkscreen scale won't rub off.*
- *Economical.*
- *Handy pocket clip.*



No. 170MM (1 - 8)



Available Models:

170MM 1 - 8 LPM
 170MM-3 1/8 - 2 1/2 LPM
 170MM-15 2 - 15 LPM

Accuracy:

170MM 1 - 8 LPM: $\pm 1/4$ LPM
 170MM-3 1/8 - 2 1/2 LPM: $\pm 1/8$ LPM
 170MM-15 2 - 15 LPM: $\pm .4$ LPM
 Length: 5 5/16"



MADA, INC.

DOMESTIC/INTERNATIONAL
 625 Washington Avenue, Carlstadt, NJ 07072
 Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

Note: The MadaMeters are hand held measuring devices that must be held as vertically as possible to insure reading.

Web Site: www.madamedical.com

r 9/10