



ASSESS[®] PEAK FLOW METER

A home monitoring device that measures peak expiratory flow, which is a valuable indicator of lung function. ASSESS[®] meets or exceeds National Asthma Education Program (NAEP) Technical Standards for Peak Flow Meters.

Performance Data

	ASSESS Standards	NAEP Standards
Accuracy	± 5%	± 10%
Reproducibility	± 5% or 10 L/min. whichever is greater	± 5% or 10 L/min. whichever is greater
Interdevice Variability	± 5%	± 5%

Ordering Information

Reorder Number	Product Description	Units/Case
710	ASSESS Standard Range	12
750	ASSESS Low Range	12
711	Adult Disposable Mouthpieces	100
712	Pediatric Disposable Mouthpieces	100

Calibration Data

	Measurement Range	Resolution
ASSESS Standard Range	60-880 L/min.	10 L/min.
ASSESS Low Range	30-390 L/min.	5 L/min.



MADA Medical
MADA Equipment
MADA International

MADA, INC.
DOMESTIC/INTERNATIONAL
625 Washington Avenue
Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509

TO ORDER CALL TOLL FREE: **800-526-6370**

See our Web Site: www.madamedical.com for a listing of our products.