

# MadaVac IV Aspirator

**NEW!**

from **MADA!**

## No. 179-II MadaVac IV

*Designed for a wide variety of  
medical suction applications.*

Small overall size  
for easy transport

Strong, balanced  
handle location  
for easy transport

Easy  
adjustment  
knob

Easy to clean,  
durable plastic  
shroud

Easy change  
hydrophobic filter

On/Off Switch

High performance  
vacuum pump

### No. 179-II Medi-Pump Specifications:

- Dimensions: L 11.5" x W 9.8" x H 9.8".
- Displacement: 15 LPM @ 0 mm Hg.
- Vacuum Range: 80-620 mm Hg @ sea level.
- Usage Rating: 30 minutes On, 30 minutes Off.
- Sound Level: 60 dba from 1 meter.
- Capacity: 800cc/Disposable Plastic Bottle.
- Service-free Life: 1500-2000 Hours.
- Weight: 2.9 Kg. (6.3 Lbs.)
- Vacuum Pump Type: Oil-less Piston.
- Cord Length 2 meters.
- Motor Type: AC, 120V, 60 Hz.
- Operating Temp: 32°-104°.

### Accessories:

- 178BE Elbow Fitting. (Each)
- 178B 800cc Suction Bottle. Single use with overflow trap. (Box 10)
- 178TB 6' Suction Tubing. (Box 10)
- UA1P712A Bacterial Filter (Each)
- 172VGTEST • Vacuum Test Gauge Assembly (Certified) (Each)

### Common Uses:

Hospitals, Nursing Homes, Homecare and Clinics.



MADA, INC.  
DOMESTIC/INTERNATIONAL  
625 Washington Ave., Carlstadt, NJ 07072  
Tel: 201-460-0454 • Fax: 201-460-3509  
To Order Call Toll Free: 800-526-6370

See our Web Site: [www.madamedical.com](http://www.madamedical.com)

r 9/10