


Mada

**No. 180
MadaMist 50
50 PSI
Compressor**



*No. 180 Air Compressor shown with
No. 182 J bracket dual set-up.*



No. 184 Gauge, Valve, Outlet

Mada's MadaMist 50 is designed to drive humidifiers, nebulizers, and mist tents.

Lightweight, compact and easy to carry.

It is designed for continuous operation and provides consistent and accurate particle size for all your nebulization needs.



MADA, INC.

DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

See our Web Site: www.madamedical.com



Ordering Information

No. 180 Air Compressor

No. 182 J Bracket, Gauge, Valve, Outlet

No. 184 Gauge, Valve, Outlet

No. 180 Specifications

- *Lightweight and compact, easy to carry. 15 1/2 lbs.*
- *Quiet.*
- *Capacity allows for efficient operation of aerosol administration devices such as humidifiers, nebulizers and mist tents.*
- *Oil-free air supply.*
- *Designed for continuous operation.*
- *Output: Maximum - 7 LPM at 50 PSIG.
Maximum - 20 LPM at 20 PSIG.*
- *Size: 9" long x 5" wide x 10" high.*
- *Outlet: Oxygen 9/16 - 18 DISS male threads.*
- *Power: 115V, 60 Hz, 1/8 Hp, 4.2 A.*
- *Includes 5 1/2 ft. power cord with hospital grade plug.*



MADA, INC.

DOMESTIC/INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

See our Web Site: www.madamedical.com

Distributed By: