

# Personal Best®

At last...  
a truly portable  
Peak Flow Meter!



## Product Specifications

Personal Best™ meets or exceeds National Asthma Education Program Technical Standards for Peak Flow Meters, based on American Thoracic Society Spirometry Waveform 24.

### Performance Data

	ASSESS Standards	NAEP Standards
Accuracy	± 5%	± 10%
Reproducibility	± 5% or 10 L/min. whichever is greater	± 5% or 10 L/min. whichever is greater
Interdevice Variability	± 5%	± 5%

### Calibration Data

	Measurement Range	Resolution
Personal Best	60-810 L/min.	10 L/min.
Low Range	50-390 L/min.	5 L/min.

### Dimensions

Closed Unit	
Length	6.5"
Width	2.0"
Height	0.8"
Weight	3.0 oz.

### Ordering Information

Reorder Number	Product Description	Units/Case
755	Personal Best Standard Range	12
756	Personal Best Low Range	12
757	Disposable Mouthpieces	100



**MADA, INC.**  
DOMESTIC/INTERNATIONAL  
625 Washington Avenue, Carlstadt, NJ 07072  
Tel: 201-460-0454 • Fax: 201-460-3509

See our Web Site: [www.madamedical.com](http://www.madamedical.com) for a listing of our products.

TO ORDER CALL TOLL FREE: **800-526-6370**

r 9/10