


Mada "R"/"M" Series Oxygen Click Style Regulators



No. R1835-8GB




No. R1833-15G




No. R1835GBHF-25



No. R1835-8GE



No. R1835-15G



No. M1935-8GB

Features:

- Compact aluminum or brass body.
- Tested to CGA standards.
- Dial-in flow control.
- Maximum inlet pressure 2,650 PSI.
- Sintered filter for safety.
- Internal reseating relief valve.
- 1 1/2" diameter gauge with protector.
- Large variety of regulators available.



MADA, INC.

DOMESTIC/INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072

Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

Web Site: www.madamedical.com

Mada Medical Oxygen Click Style Regulators

Standard Click Style Oxygen Regulators for Small Cylinders

PRODUCT NO.	FLOW RATE	INLET CONNECTION	OUTLET STYLE
R1835-4GB	1/64 - 4 LPM	CGA-870 for small cylinders	Barb Outlet
R1835-8GB	1/4 - 8 LPM	CGA-870 for small cylinders	Barb Outlet
R1835-8G	1/4 - 8 LPM	CGA-870 for small cylinders	DISS Male Outlet
R1835-15GB	1/2 - 15 LPM	CGA-870 for small cylinders	Barb Outlet
R1835-15G	1/2 - 15 LPM	CGA-870 for small cylinders	DISS Male Outlet
R1835-25GB	1/4 - 25 LPM	CGA-870 for small cylinders	Barb Outlet
R1835-25BR	1/4 - 25 LPM	CGA-870 for small cylinders	Barb Outlet, All Brass, Right Handed
R1835GBHF-25	1/4 - 25 LPM	CGA-870 for small cylinders	w/2 50 PSI DISS and 1 Barb Outlet, All Brass
R1835-8GE	1/8 - 8 LPM	CGA-870 for small cylinders	DISS Male Outlet (Extended)
R1835-15GE	1/8 - 15 LPM	CGA-870 for small cylinders	DISS Male Outlet (Extended)

Compact Click Style Oxygen Regulators

PRODUCT NO.	FLOW RATE	INLET CONNECTION	OUTLET STYLE
M1935-8B	1/2 - 8LPM	CGA-870 for small cylinders	Barb Outlet *no contents gauge
M1935-8GB	1/2 - 8 LPM	CGA-870 for small cylinders	Barb Outlet
M1935-15GB	1/2 - 15LPM	CGA-870 for small cylinders	Barb Outlet

Standard Click Style Oxygen Regulators for Large Cylinders

PRODUCT NO.	FLOW RATE	INLET CONNECTION	OUTLET STYLE
R1833-4G	1/64 - 4 LPM	CGA-540 for large cylinders	Barb Outlet
R1833-8G	1/4 - 8 LPM	CGA-540 for large cylinders	DISS Male Outlet
R1833-15G	1/2 - 15 LPM	CGA-540 for large cylinders	DISS Male Outlet
R1833GBHF-25	1/4 - 25 LPM	CGA-540 for large cylinders	Barb Outlet, Aluminum/Brass Construction

Mada Quality You Can Count On!

- Mada offers a wide range of oxygen regulators that are designed for hospital, homecare and EMS applications.
- Mada oxygen regulators offer the features you expect, including easy to read flow rates and dial-in click-style flow controller at the most competitive prices.



MADA, INC.

DOMESTIC/INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

Web Site: www.madamedical.com

Distributed By: