

# MADA O<sub>2</sub> THERAPY SHOULDER & CASE KITS



No. 1315ME



No. 1415MSE



No. 1514E



No. 1515AE



No. 1515ME



No. 1516AE

- All Mada Kits contain an oxygen cylinder, regulator, shoulder bag or case with a nasal cannula or mask and tube.
- All Mada Kits contain refillable oxygen aluminum cylinders that can be used for various medical applications.
- All Mada Kits are economically priced.
- All Mada Kits are lightweight and can be ordered filled with oxygen.
- Quality and design you can always expect from Mada.

*To order filled, remove suffix "E"*



**MADA, INC.**

DOMESTIC/ INTERNATIONAL  
625 Washington Avenue, Carlstadt, NJ 07072  
Tel: 201-460-0454 • Fax: 201-460-3509  
To Order Call Toll Free: **800-526-6370**

Web Site: [www.madamedical.com](http://www.madamedical.com)

# MADA O<sub>2</sub> THERAPY SHOULDER BAG/CASE KITS

PRODUCT NO.	CYLINDER	REGULATOR	FLOW RATE	BAG/CASE
1215ME	(M6) 1203ME	1309A	2 - 8 LPM	1209 Shoulder Bag
1313ME	(M7) 1303ME	1309A	2 - 8 LPM	1319A Carry Case
1313M-5E	(M7) 1303ME	1309A-5	1/2 - 5 LPM	1319A Carry Case
1314ME	(M7) 1303ME	1309A	2 - 8 LPM	1318A Carry Case
1314M-5E	(M7) 1303ME	1309A-5	1/2 - 5 LPM	1318A Carry Case
1315ME	(M7) 1303ME	1309A	2 - 8 LPM	1408 Shoulder Bag
1315M-5E	(M7) 1303ME	1309A-5	1/2 - 5 LPM	1408 Shoulder Bag
1415MSE	(M9) 1403ME	1309A	2 - 8 LPM	1408 Shoulder Bag
1415MS-5E	(M9) 1403ME	1309A-5	1/2 - 5 LPM	1408 Shoulder Bag
1415SE	(M9) 1403E	1443	2 - 8 LPM	1408 Shoulder Bag
1514E	(MD) 1503E	1443	2 - 8 LPM	1520 Carry Case
1514ME	(MD) 1503ME	1309A	2 - 8 LPM	1520 Carry Case
1515E	(MD) 1503E	1443	2 - 8 LPM	1509 Shoulder Bag
1515-5E	(MD) 1503E	1443-5	1/2 - 5 LPM	1509 Shoulder Bag
1515AE	(MD) 1502E	1335	2 - 8 LPM	1509 Shoulder Bag
1515A-5E	(MD) 1502E	1335-5	1/2 - 5 LPM	1509 Shoulder Bag
1515A-15E	(MD) 1502E	1335-15	2 - 15 LPM	1509 Shoulder Bag
1515ME	(MD) 1503ME	1309A	2 - 8 LPM	1509 Shoulder Bag
1515M-5E	(MD) 1503ME	1309A-5	1/2 - 5 LPM	1509 Shoulder Bag
1516AE	(MD) 1502E	M1935-8GB	1/2 - 8 LPM	1509 Shoulder Bag

Other valve combinations available upon request.

Filled Kits available upon request.



**MADA, INC.**

DOMESTIC/ INTERNATIONAL

625 Washington Avenue, Carlstadt, NJ 07072

Tel: 201-460-0454 • Fax: 201-460-3509

To Order Call Toll Free: **800-526-6370**

Web Site: [www.madamedical.com](http://www.madamedical.com)

Distributed By: