



MADA Patient Aids

Powder Coated Steel Folding Bath Bench



3310B

FEATURES:

- ▲ One piece, blow molded seat with textured surface.
- ▲ Wide base with large rubber feet for extra stability.
- ▲ 1" diameter, white powder coated steel framework.
- ▲ Built-in carry handle.
- ▲ Convenient folding mechanism for easy storage.
- ▲ No assembly required.

SPECIFICATIONS:

WEIGHT CAPACITY: 250 lbs.

SEAT WIDTH: 18³/₄".

SEAT DEPTH: 11³/₄".

HEIGHT: 17.5" from seat to floor.

PACKED: 2/cs.

