



# MADA Patient Aids

## Fold-Down Bedside Commode

**3550C**



### FEATURES:

- ▲ Heavy duty 1" diameter steel tubing.
- ▲ Folding feature provides space saving convenience for storage and transport.
- ▲ Powder coated finish.
- ▲ Easily removable back rest.
- ▲ Snap on solid plastic seat with smooth finish for easy disinfecting.
- ▲ Full length textured armrests with molded form for safe and secure grip.

### SPECIFICATIONS:

- ▲ Width between armrests 18".
  - ▲ Height adjustment from floor to hand grips 24" - 30".
  - ▲ Top of seat to floor 16½"-22½".
  - ▲ Seat Depth 15".
  - ▲ Width between legs 19¼"
- PACKING:**
- 3550C (2/cs).  
includes toilet seat and lid,  
pail and lid and splash guard.

**WEIGHT LIMIT:** 300 lbs.

### OPTIONS:

Toilet seat lid  
model 3503 (1 each).

Toilet Seat model 3504 (1 each).

Economical pail and lid  
model 3501 (12/cs).

Economical splash shield  
model 3502 (14/cs).