



# MADA Patient Aids

## Heavy Duty, Drop Arm Fold-Down 3-in-1 Commode

3555C



### FEATURES:

- ▲ No tools needed for assembly.
- ▲ Heavy duty 1" diameter steel tubing with "X" frame provides maximum stability.
- ▲ Folding feature provides space saving convenience for storage and transport.
- ▲ Powder coated - white finish.
- ▲ Easy removable back rest.
- ▲ Easy release drop arm allows for safe side transfer.
- ▲ Snap on solid plastic seat with smooth finish for easy disinfecting.
- ▲ Textured armrests with molded form for safe and secure grip.

### SPECIFICATIONS:

- ▲ Distance between armrests 22 $\frac{1}{4}$ ".
- ▲ Height adjustment from floor to hand grips 22" - 24 $\frac{3}{4}$ ".
- ▲ Top of seat to floor 16 $\frac{1}{4}$ " - 19".
- ▲ Seat depth 17".
- ▲ Width between legs 20 $\frac{3}{4}$ ".

### PACKING:

3555C (2/cs).  
Includes toilet seat with lid, pail/lid  
and splash guard.

**WEIGHT LIMIT:** 375 lbs.

### OPTIONS:

Snap on toilet seat lid (large)  
model 3505.

Toilet seat model 3504.

Economical pail and lid  
model 3501 (12/cs).

Economical splash shield  
model 3502 (14/cs).