




# MADA Patient Aids

## Patriot Two-Button Folding Walker Junior/Adult In One



### FEATURES:

1. High density vinyl hand grips.
2. High strength 1" lightweight aluminum tube construction with anodized silver finish.
3. Height adjustable 28" - 38" with push-button spring lock security.
4. Non-skid tips.
5. Two-button folding mechanism.
6. Double cross brace.
7. Optional universal support 5" wheel extensions with skis (3802PK).
  - ▲ Adjustment height 30" - 37" with 5" wheel package..
  - ▲ Nylon glide caps only (Pair) (3805GC). (Not Shown).
8. Optional nylon skis (3805NSA).

**COLOR:** Silver.

### ALSO AVAILABLE:

3825WSL-5 with 5" wheel and ski package.

3825WSL-5A with 5" wheels and skis attached

**WIDTH:** 18"

**WEIGHT LIMIT:** 300 lbs. Walking Weight.

**PACKED:** 2/cs.