



# MADA Patient Aids

## Patriot Reciprocal Folding Walker and Quad Cane



3901RWDB



"The walker  
that walks with you"



3601-7C

### FEATURES: 3901RWDB

- ▲ Padded foam hand grips.
- ▲ Height adjustable from 31" - 38" with snap pin security.
- ▲ High strength 1" lightweight anodized aluminum tube construction.
- ▲ Non-skid tips.
- ▲ Protected spring pin folding mechanism.
- ▲ Double front cross brace.

**WEIGHT LIMIT:** 275 lbs. Walking Weight.

**PACKED:** 2/cs. **COLOR:** Bronze.

### FEATURES: QUAD CANES:

- ▲ 3601-7C small base quad cane (5" x 7" base) chrome.

**WEIGHT LIMIT:**  
250 lbs. Walking Weight.

**PACKED:** 4/cs.