



MADA Patient Aids

Liberty "Excel" Lightweight Aluminum Stroll Walker



FEATURES:

1. Ultra light touch locking loop hand brakes.
2. Adjustable handle height 32.5" - 36.5" with star knob locking system.
3. Anatomically formed hand grips.
4. Extra large seat, fully padded with velcro pocket.
5. Removable shopping basket with handle and dual position hooks.
6. Casters: 8" solid tires front and rear.
7. Removable contoured padded backrest.
8. Plastic capped screws.

WEIGHT: 19 lbs.

WEIGHT LIMIT: 300 lbs.
Walking Weight.

COLOR: Burgundy.

PACKED: 1/cs.

